
<https://WKmed.by/>
+375-29-612-93-03

RIKTA

T1/T2

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DEAR CUSTOMER!

We congratulate you on purchasing the
EMITTER T1 / T2
and we express confidence that it will help
to maintain your health for many years!

The emitter is manufactured at the enterprise of JSC „MILTA – PKP GIT“, whose
quality management system meets the requirements of international standard
ISO 13485:2016.

The emitter is registered in the State Register of Medical Devices of the Russian
Federation as part of the RIKTA magneto-infrared laser therapeutic device.
Registration certificate No.ФСР 2011/10170.
Validity period — indefinitely.

CONTENTS

RU	3
EN	
1. GENERAL INFORMATION	10
2. INTENDED USE.....	10
3. SAFETY PRECAUTIONS	11
4. DELIVERY SET	12
5. DEVICE.....	12
6. TECHNICAL SPECIFICATION.....	13
7. OPERATING INSTRUCTION	13
8. STORAGE AND TRANSPORTATION.....	13
9. MANUFACTURER'S WARRANTY.....	14
FR	15
DE	21
WARRANTY CARD.....	27
WARRANTY REPAIRS.....	27

EN

1. GENERAL INFORMATION

- The present operating manual (further – the Manual) is the document that certifies basic parameters and technical specifications of the emitters of models T1-02, T1-03, T1-04, T1-05; T2-02, T2-03, T2-04, T2-05 (further emitter T1 / T2), guaranteed by the Manufacturer, and is intended for acquaintance with the operating rules of the Emitter.

The Operating manual should be kept during all the Emitter service life.

- It is necessary to read the Manual carefully before using the emitter. Particular attention should be paid to compliance with the safety precautions listed on page 5.
- The methods of carrying out treatment procedures are stated in the Methodical manual on application of RIKTA devices (further – the Methodical manual) or in individual doctor's prescription.
- The manufacturer reserves the right to make adjustments and amendments to this Manual and technical documentation at any time without notice.
- The manufacturer reserves the right to make changes in the configuration and design of the emitter, which will not affect its operation and basic technical specifications.

2. INTENDED USE

- The Emitter is intended for therapeutic noninvasive (without breaking the skin integument) effect on the patient's pathology zones.
- The Emitter should be used jointly with the magnetic-infrared laser therapeutic devices of RIKTA series. Depending on the modification, the Emitter is intended for connection only to a specific model of the magnetic-infrared laser therapeutic device of RIKTA series (the modification is specified in the nameplate on the emitter case):

T1-02 – with RIKTA -02/1 («QUANTERRA») device;

T2-02 – with RIKTA -02/1 («QUANTERRA») device;

T1-03 – with RIKTA -03/2 device;

T2-03 – with RIKTA -03/2 device;

T1-04 – with RIKTA -04/4 device;

T2-04 – with RIKTA -04/4 device;

T1-05 – with RIKTA -05 device;

T2-05 – with RIKTA -05 device.

3. SAFETY PRECAUTIONS

Using the Emitter is allowed only after acquaintance with this Manual and the Methodical manual.

- The Emitter is intended to work jointly with devices of RIKTA series (see page 10). It is not acceptable to connect the Emitter to a device for which it is not intended!
- Don't open up the Emitter and don't attempt to disassemble or modify its internal components.
- During operation, it is recommended to avoid the impact of the direct or specularly reflected radiation into the eyes.

Don't stare into the beam and don't conduct the direct observations using optical instruments.

- Before and after the treatment procedure, wipe the Emitter with gauze or cotton, moistened in a disinfectant for medical devices based on 3% hydrogen peroxide solution.

Avoid any contact of the device with alcohol-containing means! Disinfection of the Emitter should be carried out only after switching off and disconnecting the device from the power supply.

The full list of safety precautions is specified in the Operating manual for your RIKTA device. Indications and contraindications to laser therapy treatment are stated in the Methodical manual.

EN


The Emitter disposal shall meet the conventional standards and local regulations for the environment protection.

4. DELIVERY SET

Delivery set of the Emitter T1/ T2

Description	Quantity, pcs
Emitter T1 / T2	1
Operating manual	1
Package	1

5. DEVICE

EN

6 - nameplate:

6. TECHNICAL SPECIFICATION

Radiation wavelength, μm :

- laser infrared radiation..... 0,800...0,910
- broadband infrared radiation..... 0,860...0,960
- visible red radiation..... 0,600...0,700

Laser radiation pulse power

- **Emitter T1**, W..... 16^{+2}_{-4}
- **Emitter T2**, W..... 20^{+2}_{-4}

Infrared radiation average power, mW..... 60 ± 30

Red radiation average power, mW..... 7^{+3}_{-5}

Magnetic induction, mT..... 35 ± 10

Emitting window area, cm^2 $4 \pm 0,4$

Laser safety by IEC 60825-1.....class1M

Electrical safety by IEC 60601-1..... type BF

Service life of the Emitter, years.....5

7. OPERATING INSTRUCTION

 **Before connecting, make sure that the Emitter corresponds to your model of RIKTA device (see page 10).**

- Connect the cable [5] of the Emitter to the appropriate socket of the electronic control unit of your RIKTA device.
- Turn on the device.

The operating instructions are described in the Operating manual of your RIKTA device.

EN

8. STORAGE AND TRANSPORTATION

- Each Emitter is packed in the Manufacturer's package.
- The Emitter should be stored in the Manufacturer's package in the heated room at temperature in the range of $+10...+35^\circ\text{C}$ and relative humidity of no more than 80% (at $+25^\circ\text{C}$). The storage room should be free of dust, mercury vapor, acids and alkalis provoking corrosion.
- The Emitter should be transported by all types of the closed transport in Manufacturer's package. It is necessary to ensure a stable position of the package during transportation.

9. MANUFACTURER'S WARRANTY

- The Manufacturer guarantees that the Emitter parameters correspond to the values specified in the Manual, within 24 months from the date of sale of the Emitter to the consumer (if the consumer complies with the conditions of use, transportation and storage). Within the aforesaid period, the manufacturer provides free repair of the device.
- The manufacturer is not responsible for and does not make free repairs of malfunctions caused by the fault of the consumer (for example, due to careless handling, transportation, non-compliance with the rules of use, as well as in case of breaking of the seal).
- This warranty does not apply to any of the following:
 - absence of the present *Manual*;
 - absence of the *Warranty card*;
 - mechanical damage including the damage during transportation;
 - breaking of the seal.
- The warranty also becomes void if during the warranty period repair of the faulty Emitter was carried out by unauthorized persons and organizations.
- After the warranty period, the manufacturer shall repair the Emitter at the expense of the consumer.
- In the absence of the sales date and the seller's stamp in the *Warranty card*, the warranty period shall be calculated from the date of manufacturing of the Emitter.
- All issues related to repairs and maintenance should be directed to the manufacturer or to the technical centers listed on the website www.rikta.ru or in a separate warranty booklet (if available).

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